experts in muscle and bone.

Leonardo Mechanograph

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

Results Leonardo KIOSK

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Baze:

Normative values for maximum power during motor function assessment of jumping among physically active Japanese. Tsubaki A, Kubo M, Kobayashi R, Jigami H, Takahashi HE; J Musculoskelet Neuronal Interact., 9(4):263-7, 2009

Reproducibility of the jumping mechanography as a test of mechanical power output in physically competent adult and elderly subjects. Rittweger J, Schiessl H, Felsenberg D, Runge M; J Am Geriatr Soc., 52(1):128-31, 2004

Cross-sectional assessment of neuromuscular function using mechanography in women and men aged 20-85 years. Dietzel R, Gast U, Heine T, Felsenberg D, Armbrecht G.; J Musculoskelet Neuronal Interact, 13(3):274-81, 2013

Intra- and inter-rater reliability of jumping mechanography muscle function assessments. Matheson LA, Duffy S, Maroof A, Gibbons R, Duffy C, Roth J; J Musculoskelet Neuronal Interact, 13(4):480-6, 2013

Reproducibility of jumping mechanography in healthy children and adults. Veilleux LN, Rauch F; J Musculoskelet Neuronal Interact., 10(4):256-66, 2010

Assessment of musculoskeletal system in women with jumping mechanography. Dionyssiotis Y, Galanos A, Michas G, Trovas G, Lyritis GP; Int J Womens Health., 1:113-8, 2010

Normative-referenced percentile values for physical fitness among Canadians. Hoffmann MD, Colley RC, Doyon CY, Wong SL, Tomkinson GR, Lang JJ; Health Rep, 30(10):14-22, 2019

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High-frequency vibration training increases muscle power in postmenopausal women. Russo CR, Lauretani F, Bandinelli S, Bartali B, Cavazzini C, Guralnik JM, Ferrucci L; Arch Phys Med Rehabil., 84(12):1854-7, 2003

Decline of specific peak jumping power with age in master runners . Michaelis I, Kwiet A, Gast U, Boshof A, Antvorskov T, Jung T, Rittweger J, Felsenberg D; J Musculoskelet Neuronal Interact., 8(1):64-70, 2008

Jumping Mechanography Safely Evaluates Muscle Performance in Older Adults Buehring B, Valentine S, Woods A, Checovich M, Krueger D, Binkley N; ASBMR 2008

Jumping mechanography: a potential tool for sarcopenia evaluation in older individuals. Buehring B, Krueger D, Binkley N; J Clin Densitom., 13(3):283-91, 2010

A Randomized, Controlled Trial of Vitamin D Supplementation upon Musculoskeletal Health in Postmenarchal Females. Ward KA, Das G, Roberts SA, Berry JL, Adams JE, Rawer R, Mughal MZ; J Clin Endocrinol Metab; 95:4643-4651, 2010

Correlation of Different Muscle Function Tests in Older Adults. Buehring B, Fidler E, Krueger D, Binkley N; ASBMR, Annual Meeting 2011, San Diego, 2011

Age-Related Changes in Physical Function in Community-Dwelling People Aged 50-79 Years. Tsubaki A, Kubo M, Kobayashi R, Jigami H, Takahashi HE; J. Phys. Ther. Sci, 22/1:23-27, 2010

Reproducibility of jumping mechanography and traditional measures of physical and muscle function in older adults. Buehring B, Krueger D, Fidler E, Gangnon R, Heiderscheid B, Binkley N; Osteoporos Int, 26(2):819-25, 2015

Mechanography performance tests and their association with sarcopenia, falls and impairment in the activities of daily living - a pilot cross-sectional study in 293 older adults. Dietzel R, Felsenberg D, Armbrecht G; J Musculoskelet Neuronal Interact, 15(3):249-56, 2015

Greater association of peak neuromuscular performance with cortical bone geometry, bone mass and bone strength than bone density: A study in 417 older women. Belavy DL, Armbrecht G, Blenk T, Bock O, Borst H, Kocakaya E, Luhn F, Rantalainen T, Rawer R, Tomasius F, Willnecker J, Felsenberg D; Bone, 83:119-126, 2015

Effect of age and sex on jumping mechanography and other measures of muscle mass and function. Siglinsky E, Krueger D, Ward RE, Caserotti P, Strotmeyer ES, Harris TB, Binkley N, Buehring B; J Musculoskelet Neuronal Interact, 15(4):301-8, 2015

Muscle Mechanography: A Novel Method to Measure Muscle Function in Older Adults. Taani MH, Kovach CR, Buehring B; Res Gerontol Nurs, 10(1):17-24, 2017

Lower Jump Power Rather Than Muscle Mass Itself is Associated with Vertebral Fracture in Community-Dwelling Elderly Korean Women. Lee EY, Lee SJ, Kim KM, Seo DH, Lee SW, Choi HS, Kim HC, Youm Y, Kim CO, Rhee Y; *Calcif Tissue Int*, 100(6):585-594, 2017

Associations of muscle force, power, cross-sectional muscle area and bone geometry in older UK men. Zengin A, Pye SR, Cook MJ, Adams JE, Rawer R, Wu FCW, O'Neill TW, Ward KA; *J Cachexia Sarcopenia Muscle*, doi:10.1002/jcsm.12198;, 2017

Feasibility and acceptability of using jumping mechanography to detect early components of sarcopenia in community-dwelling older women. Hannam K, Hartley A, Clark EM, Aihie Sayer A, Tobias JH, Gregson CL; *J Musculoskelet Neuronal Interact*, 17(3):246-257, 2017

Elhakeem A, Hartley A, Luo Y, Goertzen AL, Hannam K, Clark EM, Leslie WD, Tobias JH: Lean mass and lower limb muscle function in relation to hip strength, geometry and fracture risk indices in community-dwelling older women.; *Osteoporos Int*, 30(1):211-220, 2019

Stolzenberg N, Felsenberg D, Belavy DL: Postural control is associated with muscle power in post-menopausal women with low bone mass.; *Osteoporos Int*, 29(10):2283-2288, 2018

Lean mass and lower limb muscle function in relation to hip strength, geometry and fracture risk indices in community-dwelling older women. Elhakeem A, Hartley A, Luo Y, Goertzen AL, Hannam K, Clark EM, Leslie WD, Tobias JH; *Osteoporos Int*, 30(1):211-220, 2019

Are jumping mechanography assessed muscle force and power, and traditional physical capability measures associated with falls in older adults? Results from the Hertfordshire Cohort Study. Parsons CM, Edwards MH, Cooper C, Dennison EM, Ward KA; J Musculoskelet Neuronal Interact, 20(2):168-175, 2020

The prevalence of sarcopenia and relationships between muscle and bone in ageing West-African Gambian men and women. Zengin A, Jarjou LM, Prentice A, Cooper C, Ebeling PR, Ward KA; J Cachexia Sarcopenia Muscle, 9(5):920-928, 2018

Age-Related Loss of Muscle Mass, Strength, and Power and Their Association With Mobility in Recreationally-Active Older Adults in the United Kingdom. Maden-Wilkinson TM, McPhee JS, Jones DA, Degens H; J Aging Phys Act, 23(3):352-60, 2015

Feasibility and Effects on Muscle Function of an Exercise Program for Older Adults. Minett MM, Binkley TL, Holm RP, Runge M, Specker BL; Med Sci Sports Exerc, 52(2):441-448, 2020

Correlation Between Mechanography and Clinical Parameters at Six Months After Hip Fracture Surgery. Min TJ, Cho J, Ha YC, Lim JY, Kang SH, Kim DK, Seo KM, Beom J; Ann Rehabil Med, 43(6):642-649, 2019

Boli:

The effect of physiotherapy in addition to testosterone replacement therapy on the efficiency of the motor system in men with hypogonadism. Baceviciene R, Valonyte L, Ceponis J.; Medicina (Kaunas), 49(2):71-7, 2013

Reduced bone strength and muscle force in women 27 years after [anorexia nervosa](#). Mueller SM, Immoos M, Anliker E, Drobnjak S, Boutellier U, Toigo M; J Clin Endocrinol Metab, :jc20151011, 2015

Oscillatory whole-body vibration improves exercise capacity and physical performance in [pulmonary arterial hypertension](#): a randomised clinical study. Gerhardt F, Dumitrescu D, Gartner C, Beccard R, Viethen T, Kramer T, Baldus S, Hellmich M, Schonau E, Rosenkranz S; Heart, 103(8):592-598, 2017

A randomized exploratory phase 2 study in patients with [chemotherapy-related peripheral neuropathy](#) evaluating whole-body vibration training as adjunct to an integrated program including massage, passive mobilization and physical exercises. Schonsteiner SS, Bauder Missbach H, Benner A, Mack S, Hamel T, Orth M, Landwehrmeyer B, Sussmuth SD, Geitner C, Mayer-Steinacker R, Riester A, Prokein A, Erhardt E, Kunecki J, Eisenschink AM, Rawer R, Dohner H, Kirchner E, Schlenk RF; Exp Hematol Oncol, 6:5, 2017

What's the secret behind the benefits of whole-body vibration training in patients with [COPD](#)? A randomized, controlled trial. Gloeckl R, Jarosch I, Bengsch U, Claus M, Schneeberger T, Andrianopoulos V, Christle JW, Hitzl W, Kenn K; Respir Med, 126:17-24, 2017

The preventive effect of sensorimotor- and vibration exercises on the onset of Oxaliplatin- or vinca-alkaloid induced peripheral neuropathies - STOP. Streckmann F, Balke M, Lehmann HC, Rustler V, Koliymitra C, Elter T, Hallek M, Leitzmann M, Steinmetz T, Heinen P, Baumann FT, Bloch W: BMC Cancer, 18(1):62, 2018

Rupp T, Butscheidt S, Jahn K, Simon MJ, Mussawy H, Oheim R, Barvencik F, Amling M, Rolvien T: Low physical performance determined by chair rising test muscle mechanography is associated with prevalent fragility fractures. Arch Osteoporos, 13(1):71, 2018

Hartley A, Gregson CL, Hannam K, Deere KC, Clark EM, Tobias JH: Sarcopenia is negatively related to high gravitational impacts achieved from day-to-day physical activity. J Gerontol A Biol Sci Med Sci, doi:10.1093/gerona/glx223:1–8, 2018

Kneis S, Wehrle A, Ilaender A, Volegova-Neher N, Gollhofer A, Bertz H: Results From a Pilot Study of Handheld Vibration: Exercise Intervention Reduces Upper-Limb Dysfunction and Fatigue in Breast Cancer Patients Undergoing Radiotherapy: VibBRa Study. Integr Cancer Ther, 17(3):717-727, 2018

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The impact of sensory and/or sensorimotor neuropathy on lower limb muscle endurance, explosive and maximal muscle strength in patients with **type 2 diabetes mellitus**. Van Eetvelde BLM, Lapauw B, Proot P, Vanden Wyngaert K, Celie B, Cambier D, Calders P; J Diabetes Complications, 107562, 2020

Whole-body vibration impedes the deterioration of postural control in patients with **multiple sclerosis**. Krause A, Lee K, Freyler K, Buhner T, Gollhofer A, Ritzmann; Mult Scler Relat Disord, 31:134-140, 2019

Jumping Joints: The Complex Relationship Between **Osteoarthritis** and Jumping Mechanography. Shere C, Fuggle NR, Edward MH, Parsons CM, Jameson KA, Cooper C, Dennison EM, Ward KA; Calcif Tissue Int, 2019

Echilibru:

Balance Performance across the Lifespan Assessed by the Leonardo Mechanograph®: A Cross-Sectional Study. Wiegmann S, Felsenberg D, Gast U, Börst H, Ambrecht G, Dietzel R; J. Funct. Morphol. Kinesiol., 5(1)

Balance Assessment in Joint **Hypermobility Syndrome**. McCarthy I; JHS, 2011

Jumping Mechanography and Classic Muscle Function Tests Effects of Age and Gender. Fidler E, Buehring B, Libber J, Checovich M, Krueger D, Binkley N; ISCD, 2012

Muscle Function, but Not Muscle Mass, is Related to Balance Confidence in **Older Adults**. Buehring B, Fidler E, Libber J, Krueger D, Binkley N; ISCD, :, 2012

Muțumesc