

HYPERTUBE

Asigurarea unui mediu curat și sănătos pentru oameni este o prioritate de top în centrul filosofiei companiei noastre, la fel cum ne-am străduit întotdeauna pentru siguranța deplină a profesioniștilor din domeniul sănătății și a pacienților noștri. Paturile de spital și cărucioarele pentru pacienți sunt articole cu risc ridicat pentru care fiecare unitate medicală trebuie să ofere îngrijirea potrivită.



HYPER TUBE asigură decontaminarea întregului pat medical, inclusiv a saltelei, într-un singur ciclu. Rapidă și fără legătură cu transportul în afara unității, soluția protejează nu numai pacienții, ci și personalul de contactul cu un pat contaminat.

Technical data

Mechanical data

Length	3150 mm
Height	2500 mm
Width	2000 mm
Housing colour	Silver bead
Safety	Motion sensor, auto off switch, app, emergency stop button
Warranty	3 years warranty (extendable)

Electrical data

Power - System	depending on project
Voltage	220 - 240 VAC; 50/60 Hz
Lamp type	Philips TUV T5
UV wavelength	253.7 nm
Power supply	Philips
Lamp life	Philips 9000 h
Additional functions	IoT + web server
Warranty	3 years warranty (extendable)

